PHOTORECEPTOR BAG

Destination of Use:

To protect the light sensetive drugs inside of bottles or mediflexes from sunlight.





TINNIT

Material:

The material is polyethilene or its kind, and also can carry up to 3 kg weight.

Dimensions:

From the cap to the end the lenght is: 330 mm (+,- 2mm) and its width is: 200mm (+,- 2 mm). The width of the cap is: 35mm.

Characteristics:

KThere is a sticker on the edge of the cap to secure the bag. There is a hole both on the bag and on its cap to hang the bag safely. There is a control line which allows to see the flow rate of the drug. This control line is also laminated to protect the drug form sunlight.

This product is disposable

| REFERENCE NUMBER | FEATURE |
|------------------|---------|
| MIP - MB - 500 | 500 сс |
| MIP - MB - 1000 | 1000 cc |

Brand : Melekırmak

